



CRYSTALYX® holder systems



CRYSTALYX® holder without a fastener clip for the 22.5 kg tub

- Can be mounted horizontally
- Allows free access to the feeding trough or feeding area
- Very durable



CRYSTALYX® holder with a fastener clip for the 22.5 kg tub

- Can be mounted horizontally or vertically
- Vertical attachment prevents a contamination of the product (e.g. through faeces)
- Locking mechanism holds the bucket firmly in place
- Allows free access to the feeding trough or feeding area
- Very durable



CRYSTALYX® holder for the 5.0 kg Mentholyx tub

- Can be mounted horizontally
- Allows free access to the feeding trough or feeding area
- The licking block is placed in the blue holder without packaging
- The blue holder can be removed and the frame can be folded, when the holder is not used.

Installation instructions:

- The holder can be mounted at the shed partition or shed wall. Install the holder close to the feeding trough so the Crystalyx® tubs can be replaced easily.
- To ensure easy access to the Crystalyx® product mount the holder according to the height of the animal.
- Crystalyx® is a natural product. When attached the holder with a fastener clip vertically and during hot temperatures leakage should not be excluded. With a daily Crystalyx® - consumption of at least 1.5 kg (200 g/day = 8 animals per pen) leakage may be avoided.